


It is important that the ear tip seals the ear canal correctly in order to achieve the attenuation level. We advise you review the diagrams and user instructions carefully, before placing in the ear.



SNR	H	M	L	Sound Technology by DYNAMIC EAR COMPANY
17 dB	20 dB	14 dB	10 dB	

Frequency (Hz)	63	125	250	500	1000	2000	4000	8000
Mean attenuation (dB)	8.8	9.0	12.2	13.2	17.1	24.2	24.8	28.0
Std. attenuation (dB)	3.8	3.1	3.7	3.1	2.5	3.5	2.3	5.4
APV (dB)	5.0	5.9	8.5	10.1	14.6	20.7	22.5	22.6

Available Ear Tip Size	Nominal Diameter
M (Medium) to L (Large)	7 - 12 mm

For additional information visit www.AheadDrumsticks.com or contact
AHEAD PRODUCTS INC. 9420 Reseda Blvd., PMB 350, Northridge, CA 91324

USER INSTRUCTIONS

WARNINGS

- If these instructions and warnings are not adhered to, the protection afforded by the product will be severely impaired.
- The product is intended to be placed in the ear and to protect the wearer from hazardous noise levels. The product is not intended for any other use.
- The product must be fitted, adjusted and maintained in accordance with the instructions. Improper fitting can reduce its effectiveness in attenuating noise.
- Performance of the product may be reduced after contact with sharp objects e.g. by puncturing or cutting.
- Regularly inspect the product to ensure serviceability.
- The product may be adversely affected when coming into contact with chemicals or alcohol.
- After first use the product can be used for up to four weeks or until soiled.
- Please ensure that the product and both hands are clean before insertion in order to reduce the risk of infection. Sharing the product can increase the risk of infection.
- Choking hazard - Keep the product out of reach of young children.
- Sudden or fast removal of the product out of the ear canal may damage the eardrum.
- Suitable hearing protection must be worn at all times in noisy surroundings.
- You are responsible for evaluating whether the use of the product affects the ability to safely perform an activity. Stop if you find it distracting or disruptive.

INSERTING INTO YOUR EARS

Roll and compress the foam part of the product between your fingers (Image A). To ease insertion, pull the top of the ear outward and upward. This is best achieved by holding the ear with your opposite hand, your arm placed around the back of the head (Image B). Gently push and twist the foam part into the ear until it sits comfortably in the ear canal. Whilst holding in position, allow the foam to expand for at least 45 seconds in order to form a good seal with the ear canal.

When the product is in position, tap gently with a finger on the outside of the hole in the plastic part. If you hear a bass drum sound effect, this is an indication of correct sealing. If placed too deeply in the ear canal, the expanding pressure of the foam may cause some discomfort. Repeat the fitting procedure if needed.

REMOVING FROM YOUR EARS

Remove the product with a slow twisting motion to gradually break the seal with the ear canal. This also eliminates any discomfort that may be felt during removal.

CLEANING AND STORAGE

Clean the product after each use to remove earwax (cerumen) and other debris. The product may be cleaned with a water dampened clean cloth. Do not clean with alcohol-based cleaning solutions, soaps or detergents. Allow to dry before use.

Always store the product in a clean resealable protective case or bag. Do not store with other objects. Recommended storage temperature range is between 14° F and 122° F.

LIMITED WARRANTY

Dynamic Ear Company (DEC) warrants its products to be free from defects in materials and workmanship for a period of two (2) years from date of purchase. If you believe this product has defects in materials or workmanship, cease use immediately and contact customer service for help.